

Buaranteed Analysis per Scoop (800 mg)
Crude Protein (Min) 5.5%
Crude Fat (Min)3.59
Crude Fiber (Max)53.59
Moisture (Max)59
Papain (Papaya)300,000 FCC-PL
Alpha-Amylase (Porcine Pancreas) 2,500 U USF
Protease (Porcine Pancreas) 2,500 U USF
Alpha-Amylase (Aspergillus oryzae) 500 DL
Lipase (Porcine Pancreas) 200 U USF
Cellulase (Trichoderma longibrachiatum) 150 CL
Hemicellulase (Aspergillus niger) 100 HCL
Xylanase (Trichoderma longibrachiatum)60 XL
Bromelain [Pineapple (Stem)] 60 GDL
Not recognized as an essential nutrient by the

AAFCO Dog or Cat Food Nutrient Profiles.

Calorie Content (Calculated): 1684 kcal ME/kg or 1.34 kcal ME/800 mg scoop.

Ingredients: Microcrystalline Cellulose, Betaine Hydrochloride, Papaya, Pineapple Stems, Porcine Trichoderma longibrachiatum, Aspergillus niger).

Pancreas, Fungal Blend (Aspergillus oryzae,

Contains a source of papain which can hydrolyze protein, 1 FCC-PU of papain activity equals the amount of papain that produces 1 mcg of tyrosine from casein under standard test conditions

²Contains a source of alpha-amylase which can hydrolyze starch. 1 U USP of amylase activity equals the amount of pancreatin that hydrolyzes 0.16 μEq of glycosidic linkages ner minute under standard test conditions.

3Contains a source of protease which can hydrolyze casein, 1 U USP of protease activity equals the amount of pancreatin that produces the equivalent 15 nmol of tyrosine from casein per minute under standard test conditions.

*Contains a source of alpha-amylase which can hydrolyze starch. 1 DU of amylase activity equals the amount of alpha-amylase that will dextrinize soluble starch with the presence of a beta-amylase at the rate of 1g/h at 30° C under standard test conditions. 5 Contains a source of lipase which can hydrolyze triglycerides. 1 U USP unit of lipase activity equals the amount of pancreatin that produces 1 µEq of acid per minute under

standard test conditions. ⁶Contains a source of cellulase which can hydrolyze beta-1, 4 glucosidic bonds. CU of cellulase activity equals the amount of activity that produces a fluidity change of 1 in 5 minutes under standard test conditions.

Contains a source of hemicellulase which can hydrolyze locust bean gum. 1 HCU of hemicellulase activity equals the amount of activity that produces a fluidity change of 1 in 5 minutes in a locust bean substrate under standard test conditions.

⁸Contains a source of xylanase which can hydrolyze xylose, 1 XU of xylanase activity equals the amount of enzyme that can liberate 1 µmol of xylose per minute under

standard test conditions. ⁹Contains a source of bromelain which can hydrolyze casein. 1 GDU of bromelain activity equals the amount that liberates 1mg of nitrogen from gelatin in 20 minutes under standard test conditions

bark & whiskers

carb starch breakdown



DIGESTIVE ENZYMES

FOR DOGS & CATS | NET WT. 4.23 OZ. (120 a)

Provides beneficial enzymes to support your pet's diet. Digestive enzymes help promote the absorption of essential nutrients

by increasing the digestibility of food. Feeding Directions: For best results, divide your pet's daily

serving amount equally between meals. Mix each portion with food. Animal Type Daily Serving Amount Scoop(s) Cats (2+ lbs.) 400 mg Toy Breed Dogs (up to 14 lbs.) Small Breed Dogs (15-29 lbs.) 1.2 g 11/2

3.2 a

This product is intended for intermittent or supplemental feeding only.

Distributed by: NHP, 125 SW 3rd Place Cape Coral, FL 33991 USA (877) 985-2696

Medium Breed Dogs (30-49 lbs.) 2 g

Large Breed Dogs (50 - 79 lbs.) 2.4 g Giant Breed Dogs (80+ lbs.)